

STATE-OF-THE-ART

X-RAY

# TOUCHSCREEN

CONTROL




**CPI**

# Advanced touchscreen technology is now a part of our

Settings & controls are right at your fingertips.

**A touch of the screen is all it takes.**



**15" Display, CPU and CD-Drive ... all-in-one.**

# Our Indico 100® family of X-Ray 'Generators-of-Choice'.

Built-in functionality helps you setup your Touchscreen.  
Optional features allow for customization.



## Anatomical Program (APR) Editor

Allows you to select, edit & store all your anatomical settings.

## Customization

Create your organization's 'look-and-feel' by adding your own company name, logo, and buttons (shape & colour).



## Integration

Another customized option that allows you to program your touchscreen to provide system integration & control.

## APR Language Selection


English, Italian, French, Spanish, Swedish & more ...  
all at your fingertips.



On-Screen Keyboard



System Utilities  
Screen



On-Screen Calendar

X-RAY  
**Generators  
of Choice**

Easy to use • Flexible • Powerful



Desktop Console (standard)



GenWare™ Service Software Interface



Upright Console (optional)

## The Indico 100® power supply

'The heart of the CPI system'



Communications & Power Industries Canada Inc.  
*communications & medical products division*  
45 River Drive / Georgetown, Ontario / Canada L7G 2J4  
Hot Line Telephone: 1-888-CPI-XRAY  
TEL: 905-877-0161 / FAX: 905-877-5327  
E-MAIL: [marketing@cmp.cpii.com](mailto:marketing@cmp.cpii.com)  
WEB: [www.cpii.com/cmp](http://www.cpii.com/cmp)

